AllergenA Basic Kit

REF

RA1000

96

REF

RA1003

192

REF

RA1004

480

REF

RA1005

960

REF

RA1100

REAGENTI DEL KIT - KIT REAGENTS - KITREAGENZIE

Reagenti / Reagents/ Reagenz	Quantità / Quantity / Menge				
	96 tests REF RA 1000	192 tests REF RA 1003	480 tests REF RA 1004	960 tests REF RA 1005	
MP	1 x 96	2 x 96	5 x 96	10 x 96	Pronto per l'uso / Ready for use / Gebrauchsfertig
CONJ E2	18 mL	2 x 18 mL	5 x 18 mL	10 x 18 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
SUBS	1 x 14 mL	2 x 14 mL	5 x 14 mL	10 x 14 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
WASH 20X	1 x 50 mL	2 x 50 mL	5 x 50 mL	10 x 50 mL	Concentrata 20X / 20X Concentrated / 20x konzentriert .
STOP	1 x 19 mL	2 x 19 mL	5 x 19 mL	10 x 19 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig

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13. LEGENDA SIMBOLI: vedere pagina 39

ENZYME IMMUNOASSAY FOR QUANTITATIVE DETERMINATION OF SPECIFIC IgE IN HUMAN SERUM

**FOR IN VITRO DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY**

1. CLINICAL APPLICATIONS

IgE takes part in occurrence and progress of allergic reactions. Binding allergens these antibodies obtain an ability to initiate release of a range of vasoactive substances from leucocytes, which define development of allergic symptoms. Definition and quantitative measurement of free allergen-specific IgE concentration have a great importance for diagnostics of allergic diseases and selection of an adequate therapy.

2. PRINCIPLE OF TEST

In the AllergenA Basic kit for detection of IgE antibodies specific to particular allergens two types of enzyme-linked immunosorbent assay are used (see Figure 1). Sandwich-type analysis is employed in microwells intended for standards. For this purpose two monoclonal antibodies against IgE with different epitopic specificity and conjugate of streptavidin with horseradish peroxidase are used. One of the antibodies is immobilized on the solid phase (interior surface of microwells) (Fig. 1 - MAT1); another one is conjugated with biotin ester (Fig. 1 - MAT2). The immobilization on the solid phase and binding of IgE molecules contained in standards happen during the first incubation (Fig. 1 - total IgE) with anti-IgE-biotin conjugate. The "capture"-variant of enzyme immunodetection is employed in wells used for samples of blood serum. One monoclonal antibody specific to IgE and immobilized on the solid phase (interior surface of microwells) (Fig. 1 -MAT1), biotinylated allergen instead of antibodies conjugated with biotin ester (at the choice of attending doctor) and conjugate of streptavidin with horseradish peroxidase are used for determination of IgE to specific allergen (Fig. 1 - specific allergen). Immobilization of different molecules present in blood serum takes place on the solid phase in microwells with samples. If there are IgE specific to chosen allergens, binding of biotinylated allergen occurs at once.

The excess of anti-IgE-biotin is removed from microwells with standards and the excess of biotinylated allergens is removed from microwells with samples of blood serum by means of washing.

During the second incubation the conjugate of streptavidin-peroxidase binds to biotinylated antibodies in microwells with standards and to the biotinylated allergen in microwells with samples tested. Removal of the excess of streptavidin-peroxidase conjugate is performed by decanting and washing of microwells.

During incubation with TMB colour of the solution changes. The colour grade is in proportion to the quantity of the bound streptavidin-peroxidase conjugate. Based on calibration curves with optical density value derived at 450 nm and 405 nm, the concentration of specific IgE in tested samples is defined after the optical density measurement of the solution in microwells.


Fig. 1. Reaction scheme

3. CONTENT OF THE KIT

It is designed for 96, 192,480 or 960 determinations.

Each plate is sufficient for an assay of 39 unknowns in duplicates or 78 unknowns in monoplicates, provided that all the strips are used simultaneously. These numbers are based on measuring the calibrators, control and blank (OD of TMB solution) in duplicates, what is recommended.

It is allowed to run the calibrators, control and blank in monoplicates, provided that the sufficient precision was proved by the user. In this case each plate is sufficient for an assay of 87 unknowns.

MP	Microplate: 12 breakable 8-wells strips each (total 96 wells) coated with anti-IgE monoclonal antibodies.
CONJ E2	Conjugate E2: 18 mL per vial of solution contains streptavidin conjugated with HRP. Preservative: Kathon 0.015 %. Ready for use.
WASH 20X	Washing Solution: 1 vial with 50 mL of surfactant in buffered saline, sufficient for preparation of 1000 mL of solution. 20X concentrated
STOP	Blocking Reagent: 1 vial with 19 mL of 1 N HCl solution. Ready to use
SUBS	TMB Solution: 1 vial with 14 mL of 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide. Ready to use

3.1 AllergenA Calibration Set

0-6 CAL	IgE Calibrators: 0.7 mL per vial of protein-based solution containing known IgE concentrations. For exact IgE concentrations, see vial labels. Preservative: Kathon 0.015 %. Ready to use.
CONTROL	IgE Control Serum: 0.7 mL per vial of protein-based solution containing known IgE concentrations. For exact range of IgE, see vial label. Preservative: Kathon 0.015 %. Ready to use.
CONJ E1	Conjugate E1: 8 mL per vial of solution containing anti-IgE monoclonal antibodies conjugated with Biotin. Preservative: Kathon 0.015 %. Ready for use.

3.2 Biotinylated Allergens

Allergen: 3 mL per vial of liquid or lyophilized solution (preservative: Kathon 0.015 %). Information about code and name of an allergen is reported on the label. One vial is sufficient for 26 assays.

3.3 Reagent Preparation

Allow all the reagents to reach a room temperature and then thoroughly stir.

MP Keep **Microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL, **CONTROL** IgE Calibrators and IgE Control Serum

Liquid Calibrators and Control Serum are ready for use.

Biotinylated Allergens: Reconstitute lyophilized Allergens with the volume of distilled or deionized water indicated on the labels of each vial.

Note: For the manual assay, please take into consideration the actual required volume when filling the amount of Conjugate, TMB Solution and Blocking Reagent into the tray.

CONJ E1 **Conjugate E1** is ready to use. Transfer Conjugate into the tray at 1.0 mL per strip.

CONJ E2 **Conjugate E2** is ready to use. Transfer Conjugate into the tray at 1.5 mL per strip.

WASH 20X Prepare required volume of **Washing Solution** by dilution the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH 20X** + 95 mL of distilled water.

Mix thoroughly, avoiding foaming. Keep the prepared wash buffer firmly closed.

SUBS **TMB Solution** is ready to use. Transfer **TMB Solution** into the tray at 1.15 mL per strip directly before dispensing. Protect TMB Solution from direct light.

STOP **Blocking Reagent** is ready to use. Transfer **Blocking Reagent** into the tray at 1.58 mL per strip before dispensing.

3.4 Sample Preparation

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kits (AllergenA Basic kit and AllergenA Calibration Set) is printed on the box label, expiry date for each component is printed on the respective label.

Upon receipt, kits should be stored at +2...+8 °C, preferably in the original kit box until the expiration date. Storage at up to +25 °C is allowed but for no more than 15 days.

Shelf life of the kits is 18 months.

Shelf life of the Allergens is in accordance with the labels.

If used for separate experiments, kits contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until expiration date;
- vials with Conjugate E1, Conjugate E2 and TMB Solution: at +2...+8 °C until expiration date;
- vials with concentrated Washing Solution and Blocking Reagent: at +2...+8 °C until expiration date;
- Washing Solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle;
- liquid (ready to use) Calibrators and IgE Control Serum: at +2...+8 °C until expiration date;
- reconstituted from lyophilized Calibrators and IgE Control Serum: at +2...+8 °C for no more than 2 months after opening;
- vials with liquid Biotinylated Allergens: firmly closed, at +2...+8 °C until expiration date.
- Biotinylated Allergens reconstituted from lyophilized preparations: at +2...+8 °C for no more than 2 months after reconstitution or at -20 °C until expiration date. Repeated freezing is allowable for no more than 4 times.

The precipitate is acceptable which vanishes after light shaking.

Note: *If used partially, kits should be utilized within 6 months after opening.*

5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- a set of calibrated variable precision pipettes with disposable tips;
- 8-channel calibrated variable precision pipette with disposable tips;
- microplate incubator / shaker (37 °C, shaking speed 400 to 800 rpm);
- manual or automatic equipment for rinsing wells;

- calibrated microplate reader (450 nm, 405 nm, 620 nm);
- graduated beakers and cylinder of appropriate volume;
- distilled or deionized water;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

6. WARNINGS AND PRECAUTIONS

In order to obtain reproducible results, the following rules must be observed:

- This kit is for *in vitro* diagnostic use only. Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only.
- Do not mix or use together reagents from different lots of the kit except TMB Solution, Blocking Reagent and Washing Solution.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Do not use TMB Solution, Blocking Reagent and Washing Solution supplied by other vendors.
- Use thoroughly clean glassware, free from contamination of metal ions or oxidating substances.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples and solutions; for this purpose, disposable tips should be used for each sample and reagent.
- Do not touch the bottom of the wells.
- IgE Calibrators should be measured in each separate assay. It is also recommended to measure each time analyte concentration in the Control.
- Reconstitute lyophilized reagents, if present, as described on the relative labels and Instructions for use. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
- If kit is used in several separate runs, it is necessary to take reagents from analyser immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator.
- Recalibration using calibration curve, obtained with kit of any other lot, is not allowed.
- TMB Solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.

- Utilise a suitable method for the correct identification of patient samples. Incorrect identification may cause a specificity losses of the system and wrong clinical results.

In order to avoid personal and environmental contamination, the following precautions must be observed:

- Use protective individual devices (ex.: disposable gloves, lab coats, etc.) while handling potentially infectious material as well as during the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be disposed of in accordance with established safety procedures.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.
- TMB Solution and Blocking Reagent (HCL 1N) should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.

As the kit contains irritant (**CONJ E1**, **CONJ E2**, **CAL**, **CONTROL** and allergens), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code 180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

7. SAMPLE COLLECTION AND STORAGE

Collect blood aseptically by venipuncture. After clotting, the serum is separated by centrifugation and transferred to an individual test tube. Do not use plasma, haemolysed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative. Store serum samples at +2...+8 °C for no more than 2 days; freeze at -20 °C or below for longer storage (but no more than 3 months). Avoid repeated freezing.

8. TEST PROCEDURE

The manual procedure is performed as described in paragraph 8.2. For internal quality control, it is advisable to use Control Serum at known concentration that allow assay performances to be monitored.

8.1 Assay Procedure for Automatic Test

We guarantee its applications on RADIM and/or SEAC/NEXT LEVEL automatic instruments.

While using a non RADIM or SEAC automatic instrument for microplate, it is under end user's responsibility to make sure that it was appropriately tested for ELISA kits.

While using for the procedure automatic analyser, refer to its relative manual.

Note: *If the kit is used on instrument in several separate experiments, it is necessary to take reagents out of analyser immediately after pipetting them into the wells of all plates, because liquid evaporates from vials. Put the reagents into refrigerator.*

8.2 Assay Procedure for Manual Test (see assay scheme, page 38)

A. Use A1 - A2 wells as "blank well".

B. Pipette 50 µL of **CAL**, **CONTROL** and 50 µL of patient's samples into the respective wells.

Note: *The number of replicates of each sample depends on the number of allergens to which specific IgE in the sample should be detected. What is more, one microwell is used for determination each allergen-specific IgE.*

C. Pipette 100 µL of **CONJ E1** into each well that contains Calibrators and Control Serum.

D. Pipette 100 µL of solution of Biotinylated Allergens into each well with tested samples.

E. Incubate for 60 minutes at 37 °C while shaking (500-800 rpm).

Note: *Total time of dispensing calibrators and samples should not exceed 15 minutes, of Allergens - 10 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

- F. Decant, then wash each well 5 times with 300 μ L of **Washing Solution** (prepared from **WASH 20X**). Make sure that after the last washing cycle the residual buffer is thoroughly aspirated from the wells.
- G. Pipette 150 μ L of **CONJ E2** into each well, except "blank well".
- H. Incubate for 30 minutes at 37 °C while shaking (500-800 rpm).
- I. Decant, then wash each well as described in paragraph F.
- J. Pipette immediately 100 μ L of **SUBS** into each well.
- Incubate for 15 minutes while shaking (500-800 rpm) in the dark at 37 °C.
- K. Pipette 100 μ L of **STOP** to all the wells and shake well for 1- 2 minutes at room temperature.
- L. Read OD at 450 nm, 405 nm and 620 nm (reference wavelength).

9. DATA PROCESSING

Plot two standard curve OD/concentration specific IgE in the calibrators (IU/mL). To plot the standard curve #1 (see Figure 2A), use the measurements results of calibrators 0...3 at a wavelength of 450 nm; to plot the standard curve #2 (see Figure 2B), use the measurements results of calibrators 0...6 at a wavelength of 405 nm.

Data processing is done by a computer-assisted analysis calculating mean OD of calibrators versus their respective specific IgE concentrations using 4PL fit.

Determine the specific IgE concentration of the sample as follows:

- Use standard curve #1 (see Figure 2A) if sample OD (at 450 nm) is less Or equal to the Cal 3;
- Use standard curve #2 (see Figure 2B) if sample OD (at 450 nm) is more than Cal 3.

Note: If reference wavelength (620 nm) is used, OD of calibrators 0...6 and specimens measured at 620 nm should be subtracted from those measured at 450 /405 nm: $OD = OD (450/405 \text{ nm}) - OD (620 \text{ nm})$. These results are used to plot the standard curves #1 and #2 and to determine specimens concentrations.

TYPICAL EXAMPLE OF STANDARD CURVE
Do not use for evaluation of real assay data!


Figure 2A. Typical example of standard curve at 450 nm


Figure 2B. Typical example of standard curve at 405 nm

9.1 Data Reliability (for OD Measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.090 ;
- concentration of IgE Control Serum must fall within the acceptability range that are shown on the vial label

Calibrator	Wavelength, nm	Average OD of Calibrators after blank subtraction
Cal 0	450	< 0.03
Cal 1	450	0.097 - 0.179
Cal 2	450	0.180 - 0.333
Cal 3	450	0.746 - 1.386
Cal 4	450	> 2.300
Cal 5	405	> 1.000
Cal 6	405	> 1.500

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

10.EXPECTED VALUES

Concentration of IgE, IU/mL	Class	Level of the specific IgE
< 0.36	0	Clinically insignificant
0.36-0.71	1	Very low
0.72-3.59	2	Low
3.60-17.99	3	Medium
18.00-49.99	4	High
50.00-100	5	Very high
> 100	6	Extremely high

11. PERFORMANCE CHARACTERISTICS OF THE ASSAY

11.1 Calibration-Traceability

Calibrators were standardized against the WHO International Standard Immunoglobulin E (IgE), human serum NIBSC code: 11/234.

11.2 Analytical Sensitivity (lower detection limit)

Analytical sensitivity of AllergenA Basic kit i.e. the lowest detectable concentration that can be distinguished from the patients' sample is 0.15 IU/mL. It was defined as mean of 10 replicates of Calibrator 0 plus 2 SD.

11.3 Specificity

No cross-reaction of monoclonal antibodies to IgE was detected with IgA, IgG, IgM and IgD.

11.4 Intra- and Inter-Assay Variation

For intra-assay CV determination, 8 serum samples were assayed in 12 replicates each. The results are shown below.

Allergens: egg white

Sample	Mean IgE concentration, IU/mL	Intra-assay	
		SD	CV %
1	0.68	0.008	1.129
2	0.85	0.005	0.57
3	1.67	0.023	1.356
4	3.15	0.027	0.85
5	12.2	0.063	0.515
6	30.5	0.368	1.207
7	65.4	0.468	0.716
8	96.8	2.427	2.508

For intra-assay quality control, the following Control Sera containing specific IgEs at known concentration are available:

REF RA9251 AllergenA Specific IgE Negative Control Serum

REF RA9249 AllergenA Specific IgE Positive Control Serum-Inhalants

REF RA9250 AllergenA Specific IgE Positive Control Serum-Foods

For inter-assay CV determination, 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 5 replicates. The results are shown below.

Sample	Mean IgE concentration, IU/mL			Inter- assay	
	Assay 1	Assay 2	Assay 3	SD	CV %
1	0.68	0.64	0.65	0.020	3.1
2	0.85	0.75	0.82	0.048	6.0
3	1.67	1.58	1.59	0.045	2.8
4	3.15	2.97	2.92	0.124	4.1
5	12.2	11.5	11.5	0.368	3.1
6	30.5	30.9	29.5	0.686	2.3
7	65.4	62.7	59.1	3.147	5.0
8	96.8	88.9	83.1	6.876	7.7

12. LIMITATION OF THE METHOD

A negative result of AllergenA Basic kit against antibiotics (penicillin G, penicillin V, cephalosporin, ampicillin and amoxicillin) doesn't exclude presence of a clinical hypersensitivity to these allergens. This result can be explained by initiation of the allergic reaction that takes place without participation of IgE antibodies or a wrong choice of blood sampling time (before the increase of specific IgE concentration in blood or after their concentration decrease).

A negative result of specific IgE assay against food-borne allergens and bites of venomous insects doesn't exclude possibility of occurrence of allergy to these allergens. If concentration of total IgE in blood serum is higher than 1000 IU/mL measured concentrations of allergen-specific IgE may be slightly lower than the real values.

It should be noted that IgE specific to certain allergens can bind other allergens with the same structure (e.g. birch pollen, nuts, meadow grass, tomato, ambrosia and melon etc.). This can happen of common antigenic determinants in allergens of different nature. These facts should be taken into account when the assay is used in clinics.

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