CE

## AllergenP Basic Kit

<b>REF</b>	PA1000	 96
<b>REF</b>	PA1003	 192
<b>REF</b>	PA1004	 480
<b>REF</b>	PA1005	 960
<b>REF</b>	PA1100	

**REAGENTI DEL KIT - KIT REAGENTS - KITREAGENZIIEN**

Reagenti / Reagents / Reagenzien	Quantità / Quantity / Menge				
	96 tests REF PA1000	192 tests REF PA1003	480 tests REF PA1004	960 tests REF PA1005	
<b>MP</b>	1 x 96	2 x 96	5 x 96	10 x 96	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>CONJ   E2</b>	1 x 13 mL	2 x 13 mL	5 x 13 mL	10 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>BUF   AS</b>	1 x 8 mL	2 x 8 mL	5 x 8 mL	10 x 8 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>WASH   20X</b>	1 x 50 mL	2 x 50 mL	5 x 50 mL	10 x 50 mL	Concentrata 20X / / 20X Concentrated / 20x konzentriert
<b>SUBS</b>	1 x 13 mL	2 x 13 mL	5 x 13 mL	10 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>STOP</b>	1 x 13 mL	2 x 13 mL	5 x 13 mL	10 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig

**AllergenP IgE Set di Calibrazione / AllergenP IgE Calibration Set REF PA1100**

*Da ordinare separatamente / To be ordered separately / Muss separat bestellt werden.*

Reagenti / Reagents / Reagenz	Quantità / Quantity / Menge	
<b>CAL   0   ...   6</b>	7 x 1,3 mL	Pronti per l'uso / Ready for use / Gebrauchsfertig
<b>CONTROL</b>	1 x 1,3 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>CONJ   E1</b>	1 x 10 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig

**Allergeni Biotinilati / Biotinylated Allergens / biotinylierte Allergene**

*Da ordinare separatamente / To be ordered separately / Muss separat bestellt werden.*

Reagenti / Reagents / Reagenz	Quantità / Quantity / Menge	
<b>Allergene / Allergen / Allergen</b>	1 x 3 mL	Liquido o Liofilizzato / Liquid or Lyophilized / flüssig oder lyophilisiert

**Sieri di Controllo / Control Sera / Kontrollseren***Da ordinare separatamente / To be ordered separately / Muss separat bestellt werden.*

<b>Reagenti / Reagents / Reagenzien</b>	<b>Quantità / Quantity / Menge</b>	<b>REF</b>
<b>Siero di Controllo Positivo-Inalanti / Positive Control Serum-Inhalants / Positivkontrolle - Inhalantien</b>	1 x 1 mL	<b>RA9249</b>
<b>Siero di Controllo Positivo-Alimenti / Positive Control Serum-Foods / Positivkontrolle - Nahrungsmittel</b>	1 x 1 mL	<b>RA9250</b>
<b>Siero di Controllo Negativo / Negative Control Serum / Negativkontrolle</b>	1 x 1 mL	<b>RA9251</b>

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**13. LEGENDA SIMBOLI:** vedere pagina 42.

# ENZYME IMMUNOASSAY FOR QUANTITATIVE DETERMINATION OF SPECIFIC IgE IN HUMAN SERUM

## FOR IN VITRO DIAGNOSTIC USE ONLY FOR PROFESSIONAL USE ONLY


### 1. CLINICAL APPLICATIONS

IgE takes part in occurrence and progress of allergic reactions. Binding allergens these antibodies obtain an ability to initiate release of a range of vasoactive substances from leucocytes, which define development of allergic symptoms. Definition and quantitative measurement of free allergen-specific IgE concentration have a great importance for diagnosis of allergic diseases and selection of an adequate therapy.

### 2. PRINCIPLE OF TEST

The **AllergenP Basic Kit** method is based on the immunoenzymatic capture system and uses the same solid phase - wells of microstrips sensitized with human anti-IgE antibodies - for all the analyses with the various allergens and for the calibration curves. The calibrators are pipetted into the wells of the first strips and the samples, each one in as many replicates as allergens to be analysed, are pipetted into the remaining wells. During the first incubation the anti-IgEs in solid phase (Fig. 1 - MAT1) capture the IgEs of the sample, both allergen-specific and non-allergen-specific (Fig. 1 – total IgE), (Fig. 1 – specific allergen).

After an initial wash to eliminate any possible interference from other immunoglobulins, (eg. allergenspecific IgGs), the anti-IgE-Biotin conjugate is added to the wells that had previously been incubated with the calibrators, and the various Allergen-Biotin conjugates are added to the dedicated wells previously incubated with the samples. During this second incubation the anti-IgE-Biotin conjugate binds to the IgE of the calibrators, captured onto the well, to form the solid-phase ||-anti-IgE : IgE : anti-IgE-Biotin sandwich (Fig. 1 – MAT2). The biotinylated allergens, on the other hand, are bound by the respective specific IgEs, captured during the first incubation, to form the solid-phase ||-anti-IgE : IgE : Allergen-Biotin immunocomplex (Fig. 1 – biotinylated allergen).

**Figure 1. Scheme of reaction**

After the wash the Streptavidin-Peroxidase conjugate is added to all the wells. This reacts both with the anti-IgE-Biotin conjugate, bound onto the wells of the calibrators, and with the Allergen-Biotin conjugate, bound onto the wells of the samples.

The last wash eliminates the non-reacted species; finally, the addition of the TMB allows detection of the Streptavidin-Peroxidase conjugate bound to the final solid-phase immunocomplex. The colour that develops is correlated to the concentration of the IgEs of the calibrators in the wells incubated with the calibrators and to the concentration of the specific IgEs in the wells incubated with the samples. The optical densities are read on a microplate reader; the sample absorbances are directly correlated to the concentration of IgE in the sample. The results are obtained by interpolation against the calibration curve and are expressed in terms of both units and classes of positivity.

### 3. CONTENT OF THE KIT (96 tests)

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-wells strips each (total 96 wells) coated with anti-IgE monoclonal antibodies.
<b>CONJ E2</b>	<b>Conjugate E2:</b> 13 mL per vial of solution containing streptavidin conjugated with HRP. Ready to use.
<b>BUF AS</b>	<b>Assay Buffer:</b> 8 mL per vial containing saline solution. Ready to use.
<b>STOP</b>	<b>Blocking Reagent:</b> 13 mL per vial of 1N HCl solution. Ready to use.
<b>SUBS</b>	<b>TMB Solution:</b> 3,3',5,5'-tetramethylbenzidine, 13 mL per vial of solution in citrate buffer containing hydrogen peroxide. Ready to use.
<b>WASH 20X</b>	<b>Washing Solution:</b> 50 mL per vial of surfactant in buffered saline, sufficient for preparation of 1000 mL of solution. <b>20X Concentrated.</b>

### 3.1 AllergenP IgE Calibration Set

<b>0-6 CAL</b>	<b>IgE Calibrators:</b> 1.3 mL per vial of protein-based solution containing known IgE concentrations. For exact IgE concentrations, see vial labels. Ready to use.
<b>CONTROL</b>	<b>IgE Control Serum:</b> 1.3 mL per vial of protein-based solution containing known IgE concentrations. For exact range of IgE, see vial label. Ready to use.
<b>CONJ E1</b>	<b>Conjugate E1:</b> 10 mL per vial of solution containing anti-IgE monoclonal antibodies conjugated with Biotin. Ready to use.

### 3.2 Biotinylated Allergens

**Allergen:** 3 mL per vial of liquid or lyophilized solution. Information about code and name of an allergen is written on the label. One vial is sufficient for 26 assays.

### 3.3 Reagent Preparation

Allow all the reagents to reach a room temperature and then thoroughly stir.

**MP** Keep microplate at room temperature (18-25°C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

#### **CAL CONTROL Calibrators and Control Serum**

Liquid Calibrators or Control Serum are ready for use.

**BUF AS** Assay Buffer is ready to use. Transfer Assay Buffer into the tray at 0.5 mL per strip.

**CONJ E1** Conjugate E1 is ready to use. Transfer Conjugate into the tray at 1.15 mL per strip.

**CONJ E2** Conjugate E2 is ready to use. Transfer Conjugate into the tray at 1.15 mL per strip.

**WASH 20X** Prepare required volume of Wash Solution by dilution the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH 20X** + 95 mL of distilled water.

Mix thoroughly, avoiding foaming. Keep the prepared wash buffer firmly closed.

**SUBS** TMB Solution is ready to use. Transfer TMB Solution into the tray at 1.15 mL per strip directly before dispensing.

Protect TMB Solution AK from direct light.

**STOP** Blocking Reagent is ready to use. Transfer Blocking Reagent into the tray at 1.15 mL per strip before dispensing.

**Allergens:** Reconstitute lyophilized allergens with the volume of distilled or deionized water indicated on the labels of each vial.

**Note:** For the manual assay, please take into consideration the actual required volume when filling the amount of Conjugate, TMB Solution and Blocking Reagent into the tray.

### 3.4 Sample Preparation

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

## 4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kits (Basic kit and Calibration Set) is printed on the box label, the expiry date for each component is printed on the respective label.

Upon receipt, kits should be stored at +2...+8 °C, preferably in the original kit box until the expiry date. Storage at up to +25 °C is allowed but for no more than 15 days.

Shelf life of the kits is 18 months.

Shelf life of the allergens is in accordance with the labels.

If used for separate experiments, kits contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until expiration date;
- vials with Conjugate E1, Conjugate E2, TMB Solution, Assay Buffer, concentrated Washing Solution and Blocking Reagent: at +2...+8 °C until expiration date;
- Washing Solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more 4 weeks, in a firmly closed bottle;
- vials with Calibrators and Control Serum: at +2...+8 °C until expiration date;
- vials with liquid allergens: firmly closed, at +2...+8 °C until expiration date.  
**Biotinylated Allergens** reconstituted from lyophilized preparations: at +2...+8 °C for no more than 2 months after reconstitution or at -20°C until expiration date. Repeated freezing is allowable for no more than 4 times.

**Note:** If used partially, kits should be utilized within 9 months after opening.

## 5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- a set of calibrated variable precision pipettes with disposable tips;
- 8-channel calibrated variable precision pipette with disposable tips;
- microplate incubator (+37 °C);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm, 405 nm, 620 nm);
- graduated beakers and cylinder of appropriate volume;

- distilled or deionized water;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

## **6. WARNINGS AND PRECAUTIONS**

**In order to obtain reproducible results, the following rules must be observed:**

- This kit is for *in vitro* diagnostic use only. Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only.
- Do not mix or use together reagents from different lots of the kit except TMB Solution, Blocking Reagent and Washing Solution.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Do not use TMB Solution, Blocking Reagent and Washing Solution supplied by other vendors.
- Use thoroughly clean glassware, free from contamination of metal ions or oxidating substances.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples and solutions; for this purpose, disposable tips should be used for each sample and reagent.
- Do not touch the bottom of the wells.
- Calibrators should be measured in each separate assay. It is also recommended to measure each time analyte concentration in the control.
- Reconstitute lyophilized reagents, if present, as described on the relative labels and Instructions for use. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
- If kit is used in several separate runs, it is necessary to take reagents from analyser immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator.
- TMB Solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
- Utilise a suitable method for the correct identification of patient samples. Incorrect identification may cause a specificity losses of the system and wrong clinical results.

**In order to avoid personal and environmental contamination, the following precautions must be observed:**

- Use protective individual devices (ex.: disposable gloves, lab coats, etc.) while handling potentially infectious material as well as during the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- TMB Solution and Blocking Reagent should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water and seek medical advice.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be disposed of in accordance with established safety procedures.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.

As the kit contains irritant (**CONJ E1** **CONJ E2** **CAL** **CONTROL** **BUF AS** and allergens), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code 180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

## 7. SAMPLE COLLECTION AND STORAGE

Collect blood aseptically by venipuncture. After clotting, the serum is separated by centrifugation and transferred to an individual test tube.

Do not use plasma, haemolysed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at 2-8°C for no more than 2 days; freeze at -20°C or below for longer storage (but no more than 3 months). Avoid repeated freezing.

## 8. TEST PROCEDURE

The manual procedure must be performed as described in paragraph 8.2.

For internal quality control, it is advisable to use Control Serum at known concentration that allow assay performances to be monitored.

### 8.1 Assay Procedure for Automatic Test

We guarantee its application on RADIM and/or SEAC/NEXT LEVEL automatic instruments.

While using a non RADIM or SEAC/NEXT LEVEL automatic instrument for microplate, it is under end user's responsibility to make sure that it was appropriately tested for ELISA kits.

While using for the procedure automatic analyser, refer to its relative manual.

**Note:** *If the kit is used on instrument in several separate experiments, it is necessary to take reagents out of analyser immediately after pipetting them into the wells of all plates, because liquid evaporates from vials. Put the reagents into refrigerator.*

### 8.2 Assay Procedure for Manual Test (see assay scheme, page 43)

A. Use **A1** well as "blank well".

B. Pipette **50 µL** of **BUF AS** into all the wells, except "blank well".

C. Pipette **50 µL** of **CAL**, **CONTROL** (in case of use of control in the assay) and **50 µL** of patient's samples into the respective wells.

**Note:** *The number of replicates of each sample depends on the number of allergens to which specific IgE in the sample should be detected. What is more, one microwell is used for the determination of each allergen-specific IgE.*

D. Incubate for **60 minutes** at **37°C**.

E. Decant, then wash each well 3 times with **300 µL** of **Washing solution** (prepared from **WASH 20X**). Make sure that after the last washing cycle the residual buffer is thoroughly aspirated from the wells.

F. Pipette **100 µL** of **CONJ E1** into each well, that contains Calibrators and Control Serum (if used), except "blank well".

G. Pipette **100 µL** of Biotinylated Allergens into each well with tested samples.

- H. Incubate for **30 minutes at 37°C**.
- I. Decant, then wash each well 3 times, as described in the paragraph E.
- J. Pipette **100 µL** of **CONJ E2** into each well, except “blank well”.
- K. Incubate for **30 minutes at 37°C**.
- L. Decant, then wash each well as described in paragraph E.
- M. Pipette immediately **100 µL** of **SUBS** into each well. Incubate for **15 minutes in the dark at room temperature**.
- N. Pipette **100 µL** of **STOP** to all the wells and shake well for 1– 2 minutes at room temperature.
- O. Read OD at **450 nm, 405 nm and 620 nm** (reference wavelength).

## **9. DATA PROCESSING**

Plot two standard curve OD/concentration specific IgE in the calibrators (IU/mL). To plot the standard curve #1 (see Figure 2A), use the measurements results of calibrators 0...3 at a wavelength of 450 nm; to plot the standard curve #2 (see Figure 2B), use the measurements results of calibrators 0...6 at a wavelength of 405 nm.

Data processing is done by a computer-assisted analysis calculating mean OD of calibrators versus their respective specific IgE concentrations using 4PL fit.

Determine the specific IgE concentration of the sample as follows:

- Use standard curve #1 (see Figure 2A) if sample OD (at 450 nm) is less or equal to the Cal 3;
- Use standard curve #2 (see Figure 2B) if sample OD (at 450 nm) is more than the Cal 3.

**Note: If reference wavelength (620 nm) is used, OD of calibrators 0...6 and specimens measured at 620 nm should be subtracted from those measured at 450/405 nm:  $OD = OD(450/405\text{ nm}) - OD(620\text{ nm})$ . These results are used to plot the standard curves #1 and #2 and to determine specimens concentrations.**

**TYPICAL EXAMPLE OF STANDARD CURVE**  
*Do not use for evaluation of real assay data!*


Figure 2A. Typical example of standard curve at 450 nm


Figure 2B. Typical example of standard curve at 405 nm

### 9.1 Data Reliability (for OD Measured at 450 nm)

The test data should meet the following criteria:

- the OD value of blank (in well A1) must be  $\leq 0.090$  (only for manual assay);
- the OD value of Cal 0 should be  $< 0.06$ .

## 10. EXPECTED VALUES

Concentration of IgE, IU/mL	Class	Level of the specific IgE
<0.36	0	Clinically insignificant
0.36-0.71	1	Very low
0.72-3.59	2	Low
3.60-17.99	3	Medium
18.00-49.99	4	High
50.00-100	5	Very high
>100	6	Extremely high

## 11. PERFORMANCE CHARACTERISTICS OF THE ASSAY

### 11.1 Calibration-Traceability

Calibrators were standardized against the WHO International Standard Immunoglobulin E (IgE), human serum NIBSC code: 11/234.

## 11.2 Analytical Sensitivity (lower detection limit)

Analytical sensitivity of **AllergenP Basic kit**, that is the lowest detectable concentration that can be distinguished in the patients' sample, is 0.15 IU/mL. It was defined as the mean of 10 replicates of Calibrator 0 plus 2 SD.

## 11.3 Specificity

No cross-reaction of monoclonal antibodies to IgE was detected with IgA, IgG, IgM and IgD.

## 11.4 Intra- and Inter-Assay Variation

For intra-assay CV determination 8 serum samples were assayed in 12 replicates each. The results are shown below.

*Allergens: egg white*

Sample	Mean IgE concentration, IU/mL	Intra-assay	
		SD	CV %
1	0.58	0.03	4.91
2	0.91	0.03	2.88
3	1.23	0.06	4.55
4	2.77	0.07	2.58
5	15.09	0.55	3.64
6	30.83	1.32	4.30
7	57.27	2.59	4.52
8	90.06	1.48	1.64

For intra-assay quality control, the following Control Sera containing specific IgEs at known concentration are available:

- REF RA9251** Negative Control Serum
- REF RA9249** Positive Control Serum-Inhalants
- REF RA9250** Positive Control Serum-Foods

For inter-assay CV determination 8 serum samples were assayed 3 times by different operators, with 1-week interval. Each specimen was assayed in 5 replicates. The results are shown below.

Sample	Mean IgE concentration, IU/mL			Inter- assay	
	Assay 1	Assay 2	Assay 3	SD	CV %
1	0.58	0.61	0.59	0.01	2.50
2	0.91	0.89	0.91	0.01	1.28
3	1.23	1.35	1.33	0.06	4.94
4	2.77	2.46	2.82	0.19	7.21
5	15.09	13.59	14.0	0.78	5.48
6	30.83	29.76	30.06	0.55	1.83
7	57.27	60.08	61.25	2.05	3.44
8	90.06	89.13	94.32	2.77	3.04

## 12. LIMITATION OF THE METHOD

A negative result of **AllergenP Basic kit** against antibiotics (penicillin G, penicillin V, cephalosporin, ampicillin and amoxicillin) doesn't exclude presence of a clinical hypersensitivity to these allergens. This result can be explained by initiation of the allergic reaction that takes place without participation of IgE antibodies or a wrong choice of blood sampling time (before the increase of specific IgE concentration in blood or after their concentration decrease).

A negative result of specific IgE assay against food-borne allergens and bites of venomous insects doesn't exclude possibility of occurrence of allergy to these allergens.

If concentration of total IgE in blood serum is higher than 1000 IU/mL measured concentrations of allergen-specific IgE may be slightly lower than the real values. It should be noted that IgE specific to certain allergens can bind other allergens with the same structure (e.g. birch pollen, nuts, meadow grass, tomato, ambrosia and melon etc.). This can happen of common antigenic determinants in allergens of different nature. These facts should be taken into account when the assay is used in clinics.

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