# AllergenP Total IgE

96

PA 1002

**REAGENTI DEL KIT - KIT REAGENTS - KITREAGENZIIEN**

<b>Reagenti / Reagents / Reagenzien</b>	<b>Quantità/ Quantity/ Menge</b>	
<b>MP</b>	1 x 96	Pronta per l'uso / Ready for use / Gebrauchsfertig
<b>0-5 CAL</b>	6 x 0,5 mL	Pronti per l'uso / Ready for use / Gebrauchsfertig
<b>CONJ E1</b>	1 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>CONJ E2</b>	1 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>BUF AS</b>	1 x 18 mL	Pronti per l'uso / Ready for use / Gebrauchsfertig
<b>DIL</b>	1 x 10 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>SUBS</b>	1 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig
<b>WASH 20X</b>	1 x 50 mL	Concentrata 20X / 20X Concentrated / 20x konzentriert
<b>STOP</b>	1 x 13 mL	Pronto per l'uso / Ready for use / Gebrauchsfertig

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**14. LEGENDA SIMBOLI:** vedere pagina 38.

# ENZYME IMMUNOASSAY FOR QUANTITATIVE DETERMINATION OF TOTAL IMMUNOGLOBULIN E (IgE) IN HUMAN SERUM.

FOR IN VITRO DIAGNOSTIC USE ONLY  
FOR PROFESSIONAL USE ONLY

## 1. CLINICAL APPLICATIONS

Normally IgE concentration in serum is very low. It increases gradually from birth to teen-age. In adults normal concentration of IgE may reach 100 IU/mL. In elderly people IgE level sometimes decreases.

IgE production is essential in anti-helminthic immunity. 15-20-fold increase in IgE concentration is observed in the case of ascariasis. But in industrialized countries detection of high IgE concentrations is mainly connected with allergic diseases. Quantitative determination of total IgE has a great prognostic value. In 75% of children born from parents with allergic diseases, serum IgE concentration is >95% from superior limit of normal range for corresponding age group. Detection of high IgE concentrations in serum by enzyme immunoassay is an important tool for differentiation between allergic diseases and other pathologies with similar clinical manifestations (such as asthma, frequent respiratory diseases, chronic rhinitis and dermatitis).

Increased concentration of total IgE in serum was also reported in patients with lymphosarcoma and Hyper-IgE syndrome.

## 2. PRINCIPLE OF TEST

The **AllergenP Total IgE** method is based on the immunoenzymatic capture system and uses the same solid phase - wells of microstrips sensitized with human anti-IgE antibodies - for all the analyses with the various allergens and for the calibration curves. The calibrators are pipetted into the wells of the first strips and the samples are pipetted into the remaining wells. During the first incubation the anti-IgEs in solid phase (Fig. 1 - MAT1) capture the IgEs of the sample.

After an initial wash to eliminate any possible interference from other immunoglobulins (eg. allergenspecific IgGs), the anti-IgE-Biotin conjugate is added. During this second incubation the anti-IgE-Biotin conjugate binds to the IgE of the calibrators, captured onto the well, to form the solid-phase || -anti-IgE : IgE : anti-IgE-Biotin sandwich (Fig. 1 - MAT2).


Figure 1. Reaction scheme

After the wash the Streptavidin-Peroxidase conjugate is added to all the wells. This reacts with the anti-IgE-Biotin conjugate, bound onto the wells of the calibrators and the samples.

The last wash eliminates the non-reacted species; finally, the addition of the TMB allows detection of the Streptavidin-Peroxidase conjugate bound to the final solid-phase immunocomplex. The intensity of the colour that develops is directly proportional to the concentration of IgE in samples or calibrators. The IgE concentration in the patient sample is read from a standard curve that is processed in each assay. The results are obtained by interpolation against the calibration curve and are expressed in terms of units.

### 3. CONTENT OF THE KIT

<b>MP</b>	<b>Microplate:</b> 1 microplate with 12 breakable 8-well strips (total 96 wells) coated with anti-IgE monoclonal antibodies.
<b>0-5 CAL</b>	<b>IgE Calibrators:</b> 6 vials with 0.5 mL each of protein-based solution containing known IgE concentrations: <b>0 - 10 - 50 - 250 - 500 - 1000 IU/mL</b> . Preservative: Kathon 0.015%. For exact IgE concentrations see vial labels. Ready to use.
<b>CONJ E1</b>	<b>Conjugate E1:</b> 1 vial with 13 mL of solution containing anti-IgE monoclonal antibodies conjugated with Biotin. Preservative: Kathon 0.015%. Ready to use.
<b>CONJ E2</b>	<b>Conjugate E2:</b> 1 vial with 13 mL of solution contains streptavidin conjugated with HRP. Preservative: Kathon 0.015%. Ready to use.
<b>BUF AS</b>	<b>Assay Buffer:</b> 1 vial with 18 mL of saline solution. Preservative: Kathon 0.015%. Ready to use.
<b>DIL</b>	<b>Sample Diluent:</b> 1 vial with 10 mL of protein-based solution. Preservative: Kathon 0.015%. Ready to use.
<b>WASH 20X</b>	<b>Washing Solution.:</b> 1 vial with 50 mL of surfactant in buffered saline, sufficient for preparation of 1000 mL of solution. <b>20X Concentrated.</b>
<b>SUBS</b>	<b>TMB Solution:</b> 1 vial with 13 mL of 3,3', 5,5' - tetramethylbenzidine, solution in citrate buffer containing hydrogen peroxide. Ready to use.
<b>STOP</b>	<b>Blocking Reagent:</b> 1 vial with 13 mL of HCl 1N solution. Ready to use.

#### 3.1 Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

**MP** Keep **Microplate** at room temperature (18-25°C) for at least 30 minutes before opening the bag. Open the bag and place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

**CAL** **Calibrators** are ready for use.

**CONJ E1 Conjugate E1** is ready to use. Transfer Conjugate into the tray at 1.15 mL per strip.

**CONJ E2 Conjugate E2** is ready to use. Transfer Conjugate into the tray at 1.15 mL per strip.

**WASH 20X** Prepare required volume of **Washing Solution** by dilution the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH 20X** + 95 mL of distilled water.

Mix thoroughly, avoiding foaming. Keep the prepared wash buffer firmly closed.

**SUBS** Protect **TMB Solution** from direct light. **TMB Solution** is ready for use. Transfer TMB Solution into the tray at 1.15 mL per strip directly before dispensing.

**STOP** **Blocking Reagent** is ready to use. Transfer Blocking Reagent into the tray at 1.15 mL per strip before dispensing.

**BUF AS** **Assay Buffer** is ready to use. Transfer Assay Buffer into the tray at 1.5 mL per strip.

### 3.2 Sample Preparation

Allow samples to reach room temperature (18-25°C). Stir samples gently in order to ensure homogeneity.

If the expected Total IgE concentration in the samples is higher than in Calibrator 5, the samples should be diluted **5-fold** or **10-fold** with **DIL** **Sample Diluent** in concordance with instructions for use. An example of manual sample dilution is as follows:

#### Sample 1 (5-fold dilution):

80 µL of **DIL** + 20 µL of serum sample

#### Sample 2 (10-fold dilution):

50 µL of **DIL** + 50 µL of the **Sample 1**.

Vortex or mix thoroughly.

## 4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; expiration date for each component is printed on the respective label.

**AllergenP Total IgE kit** should be stored at 2-8°C upon receipt, preferably in the original kit box, until the expiration date. Storage at 25°C is allowed but for no more than 15 days.

Shelf life of the kit is 18 months.

If used in several separate experiments, kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at 2-8°C until the expiration date;
- vials with Conjugate E1, Conjugate E2, Assay Buffer, Sample Diluent, concentrated Washing Solution, Blocking Reagent and TMB Solution: at 2-8°C until the expiration date;
- Washing Solution prepared for use: at room temperature (18-25°C) for no more than 5 days or at +2...+8°C for no more 4 weeks, in a firmly closed bottle;
- vials with IgE Calibrators and Control: at 2-8°C until the expiration date.

## **5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED**

- a set of calibrated variable precision pipettes with disposable tips;
- 8-channel calibrated variable precision pipette with disposable tips;
- microplate incubator (37°C);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm, 405 nm, 620 nm);
- graduated beakers and cylinder of appropriate volume;
- distilled or deionized water;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

## 6. WARNINGS AND PRECAUTIONS

**In order to obtain reproducible results, the following rules must be observed:**

- This kit is for *in vitro* diagnostic use only. Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only.
- Do not mix or use together reagents from different lots of the kit except TMB Solution, Blocking Reagent and Washing Solution.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Do not use TMB Solution, Blocking Reagent and Wash Solution supplied by other vendors.
- Use thoroughly clean glassware, free from contamination of metal ions or oxidating substances.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples and solutions; for this purpose, disposable tips should be used for each sample and reagent.
- Do not touch the bottom of the wells.
- IgE Calibrators should be measured in each separate assay. It is also recommended to measure each time analyte concentration in the Control.
- Reconstitute lyophilized reagents, if present, as described on the relative labels and Instructions for use. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
- If kit is used in several separate runs, it is necessary to take reagents from analyser immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator.
- Recalibration using calibration curve, obtained with kit of any other lot, is not allowed.
- TMB Solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
- Utilize a suitable method for the correct identification of patient samples. Incorrect identification may cause a specificity losses of the system and wrong clinical results.

**In order to avoid personal and environmental contamination, the following precautions must be observed:**

- Use protective individual devices (ex.: disposable gloves, lab coats, etc.) while handling potentially infectious material as well as during the assay.
- Do not pipette reagents by mouth.

- Do not smoke, eat, drink or apply cosmetics during the assay.
- TMB Solution and Blocking Reagent (HCL 1N) should be handled with care, as well as the reagents containing Kathon (Conjugate E1, Conjugate E2, IgE Calibrators, Sample Diluent and Assay Buffer). Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be disposed of in accordance with established safety procedures.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.

As the kit contains irritant (**CONJ E1**, **CONJ E2**, **CAL**, **BUF AS** and **DIL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

- According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code 180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

## 7. SAMPLE COLLECTION AND STORAGE

Collect blood aseptically by venipuncture. After clotting, the serum is separated by centrifugation and transferred to an individual test tube.

Do not use plasma, haemolyzed (bright red) or lipemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at 2-8°C for no more than 2 days; freeze at -20°C or below for longer storage (but no more than 3 months). Avoid repeated freezing.

## 8. TEST PROCEDURE

The manual procedure must be performed as described in paragraph 8.2.

For internal quality control, it is advisable to use Control Serum at known concentration that allow assay performances to be monitored.

### 8.1 Assay Procedure for Automatic Test

We guarantee its applications on RADIM and/or SEAC automatic instruments.

While using a non RADIM or SEAC automatic instrument for microplate, it is under end user's responsibility to make sure that it was appropriately tested for ELISA kits.

While using for the procedure automatic analyser, refer to its relative manual.

*Note: If the kit is used on instrument in several separate experiments, it is necessary to take reagents out of analyser immediately after pipetting them into the wells of all plates, because liquid evaporates from vials. Put the reagents into refrigerator.*

### 8.2 Assay Procedure for Manual Test (see assay scheme, page 39)

- A. Use **A1** well as "blank well".
- B. Pipette **150 µL** of **BUF AS** into all the wells, **except "blank well"**.
- C. Pipette **25 µL** of **CAL** and patient's samples into the respective wells.
- D. Shake gently the microplate for 10 seconds being careful not to let the content come out from the wells.
- E. Incubate for **60 minutes** at **37°C**.
- F. Decant, then wash each well **3 times** with **350 µL** of **Washing Solution** (prepared from **WASH 20X**). Make sure that after the last washing cycle the residual buffer is thoroughly aspirated from the wells.
- G. Pipette **100 µL** of **CONJ E1** into each well, **except "blank well"**.
- H. Incubate for **30 minutes** at **37°C**.
- I. Decant, then wash each well **3 times** as described in paragraph F.
- J. Pipette **100 µL** of **CONJ E2** into each well, **except "blank well"**.
- K. Incubate for **30 minutes** at **37°C**.

- L. Decant, then wash each well **3 times** as described in paragraph F.
- M. Pipette immediately **100 µL** of **SUBS** into each well. Incubate for **15 minutes in the dark at room temperature**.
- N. Pipette **100 µL** of **STOP** to all the wells and shake well for **1- 2 minutes at room temperature**.
- O. Read OD at **450 nm, 405 nm** and **620 nm** (reference wavelength).

## 9. DATA PROCESSING

Plot two standard curve OD/concentration total IgE in the calibrators (IU/mL). To plot the standard curve n° 1 (see Figure 2A), use the measurements results of calibrators 0...3 at a wavelength of 450 nm; to plot the standard curve n° 2 (see Figure 2B), use the measurements results of calibrators 0...5 at a wavelength of 405 nm.

Data processing is done by a computer-assisted analysis calculating mean OD of calibrators versus their respective total IgE concentrations using 4PL fit.

Determine the total IgE concentration of the sample as follows:

- Use standard curve n° 1 (see Figure 2A) if sample OD (at 450 nm) is less or equal to the Cal 3;
- Use standard curve n° 2 (see Figure 2B) if sample OD (at 405 nm) is more than Cal 3.

*Note: If reference wavelength (620 nm) is used, OD of calibrators 0...5 and specimens measured at 620 nm should be subtracted from those measured at 450/405 nm:  $OD = OD (450/405 nm) - OD (620 nm)$ . These results are used to plot the standard curves n° 1 and n° 2 and to determine specimens concentrations.*


**TYPICAL EXAMPLE OF STANDARD CURVE*****Do not use for evaluation of real assay data!***

Figure 2A. Typical example of standard curve at 450 nm


Figure 2B. Typical example of standard curve at 405 nm

**Data Reliability**

The data should meet the following criteria:

-average blank OD (in well A1) < 0.09;

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

**10. EXPECTED VALUES**

The spread of IgE values is extremely wide both for subjects with and without allergic disease, depending also on the type of population, i.e. geographical area, sex, and age. Therefore there is a strong overlapping of total IgE levels between atopic and non atopic populations and it is very difficult to define an upper “limit of norma”. It is recommended that each laboratory establishes its own reference range for the population of interest. The available literature provides the following informations relative to total IgE concentrations in human serum.

Mean: 29.4 IU/mL.

Age (years)	IU/mL
0-3	<10
3-4	<25
4-7	<50
7-14	<100
>=15	<150

## 11. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

## 12. PERFORMANCE CHARACTERISTICS OF THE ASSAY

### 12.1 Calibration-Traceability

**AllergenP Total IgE** kit was calibrated against the WHO International Standard Immunoglobulin E (IgE), human serum NIBSC code: 11/234.

### 12.2 Analytical Sensitivity (lower detection limit)

Analytical sensitivity **AllergenP Total IgE**, i.e. the lowest detectable concentration that can be distinguished from the patients' sample, is 3.9 IU/mL. It was defined as mean of 20 replicates of Calibrator 0 plus 2 SD.

### 12.3 Specificity:

No cross-reaction of monoclonal antibodies to IgE was detected with IgA, IgG, IgM and IgD.

### 12.4 Dilution test

Dilution tests gave values of CV% in the range 4-20 %.

### 12.5 Intra- and inter-assay variation:

For **intra-assay CV** determination, 8 serum samples were run, each in 9 replicates. The results are showed below.

Sample	Mean IgE concentration, IU/mL	Intra-assay CV	
		SD	CV, %
1	16.5	1.01	7.5
2	35.6	1.53	5.3
3	94.0	4.5	4.4
4	108	5.48	6.2
5	247	4.8	4.2
6	350	6.1	3.5
7	480	20.1	6.8
8	738	23.2	5.2

For **inter-assay CV** determination 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are showed below.

Sample	Mean IgE concentration, IU/mL			Inter-assay CV	
	Assay 1	Assay 2	Assay 3	SD	CV, %
1	16.5	15.9	16.8	0.46	2.8%
2	35.6	35.6	36	0.23	0.6%
3	94	95.9	97	1.52	1.6%
4	108	112	107	2.65	2.4%
5	247	239	260	10.60	4.3%
6	350	360	341	9.50	2.7%
7	480	430	490	32.15	6.9%
8	738	706	750	22.74	3.1%

### 13. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of **in vitro** diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

Due to the heterogeneous nature of human antibodies there might be samples that do not maintain dilution parallelism.

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