Allergen A Total IgE 2000

ENZYME IMMUNOASSAY FOR QUANTITATIVE
DETERMINATION OF TOTAL IMMUNOGLOBULIN E (IgE)
IN HUMAN SERUM

Instructions for use

(for manual test – see section 9,
for analyser “Alisei” test – see section 10)

Allergen A Total IgE 2000
(manual test)

RA1006

96

IVD

Allergen A Total IgE 2000
(analyser “Alisei” test)

RA1007

96

REF

KIT REAGENTS - KITREAGENZIEN

| Reagents / Reagenzien | Quantity/ Menge | | |
|----------------------------------|------------------------|----------------|--|
| | RA 1006 | RA 1007 | |
| MP | 1 x 96 | 1 x 96 | Ready for use / Gebrauchsfertig |
| 0-5 CAL | 7 x 0.5 mL | 9 x 0.5 mL | Ready for use / Gebrauchsfertig |
| CONJ | 1 x 18 mL | 1 x 18 mL | Ready for use / Gebrauchsfertig |
| CONTROL | 1 x 0.5 mL | 2 x 0.5 mL | Ready for use / Gebrauchsfertig |
| DIL | 1 x 3 mL | 1 x 3 mL | Ready for use / Gebrauchsfertig |
| SUBS | 1 x 14 mL | 1 x 14 mL | Ready for use / Gebrauchsfertig |
| WASH 20X | 2 x 14 mL | 1 x 50 mL | 20X Concentrated / 20x konzentriert |
| STOP | 1 x 14 mL | 1 x 50 mL | Ready for use / Gebrauchsfertig |

**ENZYME IMMUNOASSAY FOR QUANTITATIVE DETERMINATION OF TOTAL
IMMUNOGLOBULIN E (IgE) IN HUMAN SERUM**

**FOR IN VITRO DIAGNOSTIC USE ONLY
FOR PROFESSIONAL USE ONLY**

1. CLINICAL APPLICATIONS

Normally IgE concentration in serum is very low. It increases gradually from birth to adolescence. In adults, normal concentration of IgE may reach 100 IU/mL. In elderly people IgE level sometimes decreases.

IgE production is essential in anti-helminthic immunity. In the case of ascariasis, 15-20-fold increase in IgE concentration can be observed. In industrialized countries detection of high IgE concentrations is mainly connected with allergic diseases. Quantitative determination of total IgE has a great prognostic value. In 75% of children born from parents with allergic diseases, serum IgE concentration is >95% from upper limit of normal range for corresponding age group. Detection of high IgE concentrations in serum by enzyme immunoassay is an important tool for differentiation between allergic diseases and other pathologies with similar clinical manifestations (such as asthma, frequent respiratory diseases, chronic rhinites and dermatites). Also, increased concentration of total IgE in serum was reported in patients with lymphosarcoma and Hyper-IgE syndrome.

2. PRINCIPLE OF TEST

EIA-Total IgE 2000 kit is a "sandwich" type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of IgE molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. IgE molecules from the serum sample bind to both immobilized antibody and anti-IgE-peroxidase conjugate (Fig. 1). Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to the IgE level in sample. During incubation with substrate the colour is developing. The intensity of the colour is directly proportional to the concentration of IgE in specimens or calibrators. The IgE concentration in the patient sample is determined from a standard curve that is created for each assay.


Fig. 1. Reaction scheme

3. CONTENT OF THE KIT

| REF | | RA1006 | RA1007 |
|-------------------|--|--|--|
| MP | Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-IgE monoclonal antibodies | 1 pcs | 1 pcs |
| 0-6 CAL | IgE Calibrators: protein-based solution containing known IgE concentrations (0 – 2000 IU/mL). For exact IgE concentrations see vial labels | Cal 0 – 1 x 0,5 mL Cal 1 – 1 x 0,5 mL Cal 2 – 1 x 0,5 mL Cal 3 – 1 x 0,5 mL Cal 4 – 1 x 0,5 mL Cal 5 – 1 x 0,5 mL Cal 6 – 1 x 0,5 mL ready to use | Cal 0 – 1 x 0,5 mL Cal 1 – 1 x 0,5 mL Cal 2 – 2 x 0,5 mL Cal 3 – 1 x 0,5 mL Cal 4 – 1 x 0,5 mL Cal 5 – 2 x 0,5 mL Cal 6 – 1 x 0,5 mL ready to use |
| CONTROL | IgE Control: protein-based solution containing known IgE concentration. For exact range of IgE concentration, see vial labels | 1 x 0,5 mL ready to use | 2 x 0,5 mL ready to use |
| CONJ | Conjugate: solution containing anti-IgE monoclonal antibodies conjugated with HRP | 1 x 18 mL ready to use | 1 x 18 mL ready to use |
| DIL | Sample Diluent: protein-based solution | 1 x 3 mL ready to use | 1 x 3 mL ready to use |
| WASH P 20X | Wash Solution P, 20X concentrated: surfactant in buffered saline | 2 x 14 mL for preparation of 2 x 280 mL of solution | 1 x 50 mL for preparation of 1000 mL of solution |
| SUBS | TMB Solution: 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide | 1 x 14 mL ready to use | 1 x 14 mL ready to use |
| STOP | Blocking Reagent: 1 N HCl solution | 1 x 14 mL ready to use | 1 x 50 mL ready to use |

Note: Extra vials of **CAL** 2, 5 and **CONTROL** in RA 1007 are provided for recalibration on reference calibration curve. For detailed information see instructions for use for automatic analyser.

4. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly mix.

[MP] Keep **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal firmly.

CAL **CONTROL**

Calibrators and control are ready to use.

CONJ Conjugate is ready to use.

[WASH P] 20X Prepare required volume of **Wash Solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

14 mL of **[WASH P] 20X** + 266 mL of water.

Mix thoroughly, avoid foaming.

[SUBS] Substarte is ready to use. Protect **TMB Solution** from direct light.

[STOP] Blocking Reagent is ready to use.

[TRIAL] Working solution of **Trial** for washing internal hydraulic lines and needles of automatic analyser “Alisei” (see section 10). Prepare required volume of **Trial Solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL of **[TRIAL] 5000X** + 9998 mL of water.

Mix thoroughly, avoid foaming.

5. SAMPLE PREPARATION

Allow samples to reach room temperature (+18...+25 °C). Mix samples gently in order to ensure homogeneity.

If expected IgE concentration in the samples is higher than in calibrator 5 (or in calibrator 6), the samples should be diluted **20-fold** with **DIL Sample Diluent** manually according to instructions for use. The example of manual sample dilution as follows:

285 µL of **DIL Sample Diluent** + **15 µL** of serum sample, vortex or mix thoroughly.

6. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; expiration date for each component is printed on the respective label.

EIA-Total IgE 2000 should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 15 days, including at +37 °C for no more than 1 day.

Shelf life of the kit is 18 months.

After initial opening, the kit is stable for 6 months if stored at +2...+8 °C.

If used in several separate experiments, kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with Calibrators, Control: at +2...+8 °C for no more than 6 months;
- vials with Conjugate, Sample Diluent, concentrated Wash Solution, Blocking Reagent and concentrated Trial Solution: at +2...+8 °C until the expiration date;
- vial with TMB Solution: at +2...+8 °C for no more than 6 months, protected from light;
- Wash Solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8 °C for no more than 4 weeks, in a firmly closed bottle;
- Trial Solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

7. WARNINGS AND PRECAUTIONS

To obtain reproducible results, follow the rules below:

- This kit is for *in vitro* diagnostic use only. Operator should follow the instruction for use closely in order to ensure reliable data.
- Use thoroughly cleaned glassware, free from contamination of metal ions or oxidizing substances.
- Use deionized or distilled water stored in clean containers.
- Use new disposable tips for each sample and reagent to avoid any contamination among samples and solutions.
- Do not mix or use together reagents from different lots of the kit except TMB Solution, Blocking Reagent and Wash Solution.
- Do not use TMB Solution, Blocking Reagent and Wash Solution supplied by other vendors.
- Do not touch the bottom of the wells.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Reconstitute lyophilized reagents, if present, as described on the corresponding labels and Instructions for use. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
- Calibrators must be used for each new assay. It is also recommended to measure each time analyte concentration in the control.
- Recalibration using calibration curve, obtained with kit of any other lot, is not allowed.
- TMB Solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
- Use a suitable method for the correct identification of patient samples. Incorrect identification may cause loss of specificity of analysis and wrong clinical results.

To avoid any personal or environmental hazards, follow the recommendations below:

- Use individual protective items (lab coats, gloves, goggles, etc.) while handling potentially infectious material and during the assay procedure.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- TMB Solution and Blocking Reagent (1 N HCL) should be handled with care, as well as the reagents containing Kathon (Conjugate, Calibrators, Control and Sample Diluent). Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can

guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be disposed of in accordance with established safety procedures.

- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.

As the kit contains irritant (CONJ DIL CAL CONTROL), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

- According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code 180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

8. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipemic (milky) serum samples as well as samples containing sodium azide as preservative. Store serum samples at +2...+8 °C for no more than 5 days. Aliquot and freeze samples for longer storage (-20 °C and lower). Avoid repeated freezing.

9. MANUAL TEST

9.1. Equipment and materials required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator/shaker (+37 °C, shaking speed 400–800 rpm);
- manual or automatic equipment for rinsing wells;
- microplate calibrated reader (450 nm, 405 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents using 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

9.2. Test Procedure (See assay scheme, page 38)

For internal quality control, it is advisable to use Control Serum at known concentration that allow assay performances to be monitored.

It is possible to do manual test in 2 options:

Option 1: If expected IgE concentration in the samples is lower than in calibrator 5, it is recommended to use CAL set: CAL 0, CAL 1, CAL 2, CAL 3, CAL 4, CAL 5. If obtained IgE concentration in sample is higher than in CAL 5, the sample should be diluted (**see section 5**) and retested.

Option 2: If expected IgE concentration in the samples is higher than in calibrator 5, it is recommended to use CAL set: CAL 0, CAL 1, CAL 2, CAL 4, CAL 5, CAL 6. If obtained IgE concentration in sample is higher than in CAL 6, the sample should be diluted (**see section 5**) and retested.

Note: It is allowed to use all CAL set (CAL 0 – CAL 6).

A. Pipette: 150 µL of Conjugate CONJ into each well except wells A1-A2 (blank).

B. Pipette: 20 µL of Calibrators CAL, Control CONTROL and patients' samples in duplicates into the respective wells.

Leave wells A1-A2 empty for blank.

Note: total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.

C Incubate strips for 60 minutes at +37 °C while shaking (400-800 rpm) or for 120 minutes at room temperature without shaking (previously pre-shake for 1-2 minutes at room temperature).

D. Decant, then wash each well 5 times with 300 µL of Wash Solution (prepared from WASH P 20X). Make sure that after the last washing cycle the residual buffer is thoroughly aspirated from the wells. It is advisable to use an automatic microplate washer.

E. Pipette 100 µL of TMB Solution SUBS into each well (including blank); incubate strips at room temperature (+18...+25 °C) in the dark for 15-30 minutes, depending on the colour intensity, or 10 minutes while shaking (400-800 rpm) at +37 °C.

F. Pipette 100 µL of Blocking Reagent STOP into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB Solution. Shake for 1-2 min at room temperature.

G. Read OD at 450 nm (and at 405 nm, if it is used option 2) within 20 min.

9.3. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = 2.28 - 0.06 = 2.22

Data processing is done by a computer-assisted analysis calculating mean OD of calibrators versus their respective total IgE concentrations using 4PL fit.

- For **Option 1** use ODs of CAL 0, CAL 1, CAL 2, CAL 3, CAL 4, CAL 5 **at 450 nm** to draw calibration plot (see typical standard curve, fig. 2).

Determine the total IgE concentration of the sample.

- For **Option 2** use ODs of CAL 0, CAL 1, CAL 2, CAL 4, CAL 5, CAL 6 **at 450 nm** to draw calibration plot (see typical standard curve, fig. 3A) and ODs of CAL 0, CAL 1, CAL 2, CAL 4, CAL 5, CAL 6 **at 405 nm** to draw calibration plot (see typical standard curve, fig. 3B).

Determine the total IgE concentration of the sample using:

- curve 3A, if OD of the sample at 450 nm is lower than OD of CAL 5.
- curve 3B, if OD of the sample at 450 nm is higher than OD of CAL 5.

Do not extrapolate IgE concentration in the sample if it is above the nominal value of CAL 5 (at 450 nm) or CAL 6 (at 405 nm) on standard curve. In this case the sample should be diluted 20-fold with **Sample Diluent** and retested. Multiply the measured concentration of pre-diluted samples by dilution factor (20-fold).

If IgE concentration in 20-fold diluted sample still exceeds **CAL 6** concentration proceed with next dilution or assign concentration "above 40000 IU/mL" to this sample.


Fig. 2. Example of typical standard curve for Option 1.
Do not use for evaluation of real assay data.


Fig. 3A. Example of typical standard curve for Option 2 (at 450 nm).
Do not use for evaluation of real assay data.


Fig. 3B. Example of typical standard curve for Option 2 (at 405 nm).
Do not use for evaluation of real assay data.

9.4. Data Reliability

The data should meet the following criteria:

- average blank OD **at 450 nm** (in wells A1-A2) ≤ 0.09 ;
- average OD of Cal 5 **at 450 nm** ≥ 1.0 (after blank subtraction);
- average OD of Cal 6 **at 405 nm** ≥ 0.5 (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown on the vial label.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

10. AUTOMATIC TEST

10.1. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “ALISEI Q.S.”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL¹;
- **TRIAL 5000X** Trial Solution, 5000X concentrated²;
- deionized or distilled water;
- latex or plastic gloves.

10.2. Test Procedure

We guarantee our kit applications on Seac/Next Level automatic instruments.

While using a non Seac/Next Level automatic instrument for microplates, it is under end user's responsibility to make sure that it was appropriately tested for ELISA kits.

While using automatic analyser for the procedure, refer to analyser's manual.

Note: If the kit is not used at once, it is necessary to take reagents out of analyzer immediately after pipetting them into the wells of all plates, because liquid evaporates from vials. Put the reagents into refrigerator.

10.2.1. Assay Procedure for Automatic Test

While using for the procedure analyser “Alisei”, refer to its corresponding manual.

The analysis on EIA analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of Total IgE 2000 determination is stored in analyser memory.

11. EXPECTED VALUES

Each laboratory should establish its own reference range of IgE concentrations. The following reference range is recommended (see below):

| Age (years) | Mean (IU/mL) |
|-------------|--------------|
| <1 | 0 - 15 |
| 1 - 6 | 0 - 60 |
| 6 - 10 | 0 - 90 |
| 10 - 16 | 0 - 200 |
| adults | 0 - 100 |

¹ Delivered by separate order

² Delivered by separate order

12. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

13. PERFORMANCE CHARACTERISTICS OF THE ASSAY

13.1 Calibration-Traceability

EIA-Total IgE 2000 kit was calibrated against the WHO 3rd International Standard for human serum IgE, NIBSC code: 11/234.

13.2 Analytical Sensitivity

Analytical sensitivity of **EIA-Total IgE 2000** kit, i.e. concentration that can be distinguished from zero calibrator, is 8 IU/mL. It is defined as mean OD of 10 replicates of Calibrator 0 plus 2 SD.

13.3 Specificity

According to Product specification, provided by the supplier, no cross-reactivity was detected between anti-IgE monoclonal antibodies used in the assay and IgG, IgM and IgA.

13.4 Measurement Range

EIA-Total IgE 2000 kit was validated for measurement of IgE concentration within the concentration range (without dilution) of 8 - 2000 IU/mL.

13.5 Measurement Units

In **EIA-Total IgE 2000** kit the concentrations of calibrators are specified in IU/mL. To convert into ng/mL, multiply the concentration in IU/mL by 2.4.

13.6 Hook Effect

For **EIA-Total IgE 2000** kit **high dose hook effect** was not detected for concentrations up to 10 000 IU/mL. **High dose hook effect** was determined by spiking Calibrator 0 with antigen.

13.7 Intra- and Inter-Assay Variation (Precision)

Intra-assay CV was determined on 8 serum samples, each in 9 replicates. The results are shown below.

| Sample | Mean IgE concentration, IU/mL | Intra-assay CV | |
|--------|-------------------------------|----------------|-------|
| | | SD | CV, % |
| HS 1 | 24,3 | 1,84 | 7,6 |
| HS 2 | 27,1 | 1,27 | 4,7 |
| HS 3 | 57,9 | 1,56 | 2,7 |
| HS 4 | 188,2 | 3,11 | 1,7 |
| HS 5 | 167,4 | 3,39 | 2,0 |
| HS 6 | 736,2 | 22,91 | 3,1 |
| HS 7 | 1096,8 | 37,90 | 3,5 |
| HS 8 | 1958,4 | 40,16 | 2,1 |

Inter-assay CV was determined on 8 serum samples, each in 9 replicates, that were assayed 3 times by different operators with 1-week interval. The results are shown below.

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14. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. A qualified physician must consider all available clinical and laboratory findings to make proper diagnosis. Due to the heterogeneous nature of human antibodies there might be samples that do not maintain dilution parallelism.

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