MBL

CE

MESACUP ECP TEST

Cat. No. 7618E: 96 wells
Cat.No. 7618E: 96 puits
Art.-Nr. 7618E: 96 Auftragsstellen
Cat. No. 7618E: 96 pocillos
Cat. n. 7618E: 96 pozzetti
Cat. N° 7618E: 96 poços
Kat. Nr. 7618E: 96 brunnar
Ap. kat. 7618E: 96 βυθισμάτων
Kat. No. 7618E: 96 kuyucuk

FOR PROFESSIONAL USE ONLY
UNIQUEMENT À USAGE PROFESSIONNEL
NUR ZUM PROFESSIONELLEN GEBRAUCH
SOLO PARA USO PROFESIONAL
SOLO PER USO PROFESSIONALE
APENAS PARA USO PROFISSIONAL
ENDAST FÖR PROFESSIONELLT BRUK
МОНО ΓΙΑ ΕΠΑΓΓΕΛΜΑΤΙΚΗ ΧΡΗΣΗ
SADECE PROFESYONEL AMAÇLI KULLANIM İÇİNDİR

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

KDX Nagoya Sakae Bldg. 10F

4-5-3 Sakae, Naka-Ku, Nagoya, Aichi, 460-0008 Japan

Tel: +81 52-238-1901 Fax : +81 52-238-1440 URL <http://www.mbl.co.jp>

- CONTENTS -
- CONTENU -
- INHALTSVERZEICHNIS -
- CONTENIDO -
- CONTENUTI -
- CONTENDO -
- INNEHÅLL -
- ΠΕΡΙΕΧΟΜΕΝΑ –
- İÇİNDEKİLER –

English	1
Français	8
Deutsch	15
Español	22
Italiano	29
Português	
Svenska	
Ελληνικά	
Türkçe	36

Symbols/Symboles/Symbole/Símbolos/Simboli/ Simbolos/Symboler/Σύμβολα/Semboller	43
---	----

English

INTENDED USE

The MESACUP ECP TEST is quantitative assay kit for human ECP in serum by sandwich ELISA.
The MESACUP ECP TEST is intended for in vitro diagnostic use.

SUMMARY AND EXPLANATION

Eosinophil Cationic protein (ECP), Eosinophil derived Neurotoxin (EDN) and Major Basic Protein (MBP) are known to be major protein-mediators derived from activated eosinophils. ECP and EDN are found in matrix of granules in eosinophils whereas MBP is found in core of granules. ECP and EDN are members of the ribonuclease A superfamily. ECP and MBP have high cytotoxicity. Activated eosinophile play an important role in the late asthmatic response and in the asthmatic airway inflammation. As ECP is secreted from activated eosinophils, ECP can be a marker of eosinophile activation and degranulation.

The MESACUP ECP TEST is for measuring ECP specifically with high sensitivity by ELISA.

PRINCIPLE

The MESACUP ECP TEST measures human ECP by sandwich ELISA.

In the wells coated with anti-human ECP monoclonal antibody, samples to be measured or standards are incubated. After washing, a peroxidase conjugated anti-human ECP polyclonal antibody is added into the microwells and incubated. After another washing, the peroxidase substrate is added into microwells and incubated for an additional period of time. An acid solution is then added to each well to terminate the enzyme reaction and to stabilize the developed color. The optical density (O.D.) of each well is then measured at 450 nm using a microplate reader. The concentration of ECP is calibrated from a standard curve based on reference standards.

BRIEF ASSAY PROCEDURE

<Sample Incubation>	Add 100 µL of diluted sample (1:5) to each well of the microwell plate.
(20~25°C) 60 min.	↓
<Wash>	Wash
	↓
<Conjugate Incubation>	Add 100 µL of Conjugate Reagent to each well.
(20~25°C) 60 min.	↓
<Wash>	Wash
	↓
<Substrate Incubation>	Add 100 µL of Substrate to each well.
(20~25°C) 10 min.	↓
<Stop Reaction>	Add 100 µL of Stop Solution to each well.
	↓
	Read Absorbance at 450 nm
	↓
	Interpretation of Result

REAGENTS AND STORAGE

Materials	Quantity (96 wells)
Microwell Strips coated with anti-human ECP antibody	8-well x 12 strips
ECP Standard 0.09% sodium azide, 1% Goat Serum Ready-to-use	1.3 mL x 6 vials
Assay Diluent 0.09% sodium azide, 1% Goat serum Ready-to-use	50 mL x 1 vial
Conjugate Reagent (Horseradish Peroxidase conjugated anti-human ECP polyclonal antibody) 1% BSA Ready-to-use	15 mL x 1 vial
Wash Concentrate (20 x concentrated solution)	100 mL x 1 vial
Substrate TMB/H ₂ O ₂ Ready-to-use	20 mL x 1 vial
Stop Solution 0.5 mol/L H ₂ SO ₄ Ready-to-use	20 mL x 1 vial
Positive Control 0.09% sodium azide, 1% Goat Serum	1.3 mL x 1 vial
Negative Control 0.09% sodium azide, 1% Goat Serum	1.3 mL x 1 vial

Both unopened and first opened kit components are stable until the labeled expiration date at 2-8°C.

PRECAUTIONS

- (1) This product is for in vitro diagnostic use only. Do not use in human beings.
- (2) Do not use the kit components beyond the stated expiration dates.
- (3) Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- (4) Avoid contact of the reagents with eyes, skin and clothing. Reagents on skin must be washed away with plenty of water. TMB is irritating and the Stop Solution consists of 0.5 mol/L sulfuric acid, which is a poison and is corrosive.
- (5) ECP Standard, Assay Diluent, Positive Control and Negative Control contain sodium azide (0.09%) as a preservative and must be handled with caution. Do not ingest or allow contact with skin or mucous membranes. Sodium azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush with plenty of water when disposing materials containing sodium azide into a drain.
- (6) ECP Standard, Assay Diluent, Conjugate Reagent, Positive Control and Negative Control contain materials of animal origin, taken from non-infected animals. These components, however, should be treated as potential biohazards in use and for disposal.
- (7) Do not mix reagents from other kits.
- (8) All reagents must be brought to room temperature (20-25°C) before starting the assay.

- (9) Do not expose the kit to direct sun during the test and storage.
- (10) Avoid microbial and cross contamination of reagents or samples.

Contamination of the TMB substrate solution with conjugate or other oxidants will cause the solution to change color prematurely.
- (11) Incubation temperatures above or below normal room temperature (20-25°C), shorter or longer incubation times and inaccurate dilutions may give erroneous results.
- (12) The wells must be properly rinsed with Wash solution to avoid false positive results.
- (13) Carefully pipette each sample and reagent to avoid cross contamination between microwells, avoid foaming.
- (14) All unused Microwell Strips should be returned to the ziplock pouch, which must be carefully resealed to avoid moisture absorption.

Do not allow wells to become dry during an assay procedure.
- (15) Wash Concentrate may become turbid at 2-8°C, this does not cause inconsistent results.

Dissolve the Wash Concentrate completely when preparing the Wash Solution.
- (16) Materials used for the test should be disposed or treated as shown below.

Soak in 2% final conc. glutaraldehyde solution for more than 1 hour, soak in 0.1% sodium hypochlorite solution (available chlorine: approx. 1,000 ppm) for more than 1 hour or autoclave at 121°C for more than 20 minutes.
- (17) The results obtained from this assay are aids in the diagnosis only. Each physician must interpret these results in light of the patient's history, physical findings, and other diagnostic procedure.
- (18) In the event of damage to the protective packaging, the kit components should be treated or disposed of in accordance with the relevant instructions and regulations of each country.

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader (wavelength: 450 nm, 620 nm/reference)
- Multichannel micropipette (e.g. 100 µL - 300 µL)
- Single channel pipette (e.g. 10 µL & 100 µL)
- Autowasher or washing bottle
- Deionized or distilled water
- One liter graduated cylinder for preparation of Wash Solution
- Test tubes for patient sample dilutions (e.g. 1,000 µL)
- Disposable pipette tips
- Paper towels
- Microplate cover

PROCEDURE

■ PREPARATION OF REAGENTS

- 1) Bring all assay materials to room temperature (20-25 °C) prior to use.
- 2) Microwell Strips: Remove required Microwell Strips from the pouch and place them in the frame. Promptly return unused strips to refrigerated storage.
- 3) Wash Solution: The Wash Concentrate must be diluted prior to use. Dilute the Wash Concentrate 1: 10 by adding 100 mL of the concentrate to 900 mL of distilled water. The diluted Wash Solution is

stable for 2 weeks at 2-8 °C

4) Do not dilute ECP Standard, Assay Diluent, Conjugate Reagent, Substrate and Stop Solution. These reagents are ready- to-use.

■ PREPARATION OF CONTROLS

Dilution of controls

Dilute Positive Control and Negative Control 1:5 with Assay Diluent.

For example, by adding 50 µL of control to 200 µL of Assay Diluent.

* Dilute controls for each assay because they cannot be stored.

■ PREPARATION OF SAMPLES

1) Sampling

Serum levels of ECP can be varied by *in vitro* effect after sampling. There is a possibility that *in vitro* release of ECP from eosinophils can occur during coagulation. Therefore, each laboratory is recommended to establish its own sampling procedure. The following sampling procedure has been reported.

Blood for measurement of ECP was collected in Vacutainer SST tube (Becton Dickinson Vacutainer Systems, Franklin Lakes, N.J., USA). The blood allowed to clot at room temperature for exactly 60 minutes. Thereafter, the serum was separated by centrifugation for 10 minutes at 1200 x g, and then transferred to a new test tube and stored at -20 °C until the analyses were performed.

2) Dilution of sample (Human serum)

Dilute each sample 1:5 with Assay Diluent.

For example, by adding 50 µL of sample to 200 µL of Assay Diluent.

* Dilute samples for each assay because they cannot be stored.

3) Storage

Use fresh patient sera. Samples should be tested as soon as possible after collection.

If storage is needed, samples should be stored at -20°C or lower. Do not repeat freezing and thawing.

■ ASSAY PROCEDURE

STEP 1. (SAMPLE INCUBATION)

Transfer 100 µL of ECP Standard, each diluted sample, Positive Control and Negative Control, with a multichannel micropipette, into the appropriate microwells. (Do not dilute ECP Standard.)

* Incubation starts on pipetting into the microwells. Pipetting should be completed as quickly as possible.

Cover the wells with a microplate cover and incubate for 60 minutes at room temperature (20-25°C).

STEP 2. (WASHING)

Aspirate or discard the well contents. Fill the well with Wash Solution and then completely aspirate or discard the contents. Wash 4 times. Tap the plate on a paper towel to remove any remaining Wash Solution. When autowasher is used, wash 4 times.

* Each laboratory is recommended to confirm its own appropriate washing times and other conditions.

* Wash Solution should be used at 20-25°C.

STEP 3. (CONJUGATE INCUBATION)

Add 100 µL of the Conjugate Reagent to each well with a multichannel micropipette.
Cover the wells with the microplate cover and incubate for 60 minutes at room temperature (20-25°C).

* Do not return the Conjugate Reagent once taken out of vial.

STEP 4. (WASHING)

Wash the microplate as described in STEP 2.

STEP 5. (SUBSTRATE INCUBATION)

Add 100 µL of the Substrate to each well with a multichannel micropipette.

* Do not return the Substrate once taken out of vial.

Cover the wells with the microplate cover and incubate for 10 minutes at room temperature (20-25°C).

STEP 6. (STOPPING REACTION)

Add 100 µL of the Stop Solution to each well with a multichannel micropipette.

■ READING

Read the absorbance of each well at 450 nm. If a dual wavelength plate reader is used, set the test wavelength at 450 nm and the reference at 620 nm.

* Reading should be done within 20 minutes after stopping the reaction.

* Before reading the plate, ensure that the bottom of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells.

■ CALCULATION OF RESULTS


Calculate the mean absorbance value of each standard and plot against log standard concentration on suitable graph paper. The concentrations of the samples can then be read from the standard curve.

If using a program to automatically calculate the concentration, it is recommended that the best fitting curve is used.

Report the ECP concentration of samples by multiplying the value read from the standard curve by dilution factor (e.g. x 5).

*An international reference material for ECP is not available. ECP was purified and measured its protein concentration by Ultraviolet absorbance method and protein measurement method.

■ EXAMPLE OF STANDARD CURVE

■ QUALITY CONTROL

Each assay result should meet the following criteria.

A₄₅₀ of ECP Standard 1: ≤ 0.2

A₄₅₀ of ECP Standard 6: ≥ 1.2

The Positive and Negative Controls must give the following results:

ECP Value (ng/mL): As described in each label.

If any of these are not met, the results are invalid and the test should be repeated.

Before repeating assay, check the following procedure.

- Incubation Temperature
- Incubation Time
- Washing

■ CONCENTRATION OF HUMAN ECP IN NORMAL HUMAN SERUM.

Serum sample from healthy blood donors were assayed by the MESACUP ECP TEST.

148 healthy donor sera were measured. After removing those samples measuring over the mean + 3SD, the new mean + 3SD was determined to be: mean + 3SD = 15.6 ng/mL.

Sample n = 148 mean = 4.06 ng/mL SD = 3.86 mean+3SD = 15.6 ng/mL

■ LIMITATIONS

As with other diagnostic test procedures, the results obtained with the MESACUP ECP TEST serve only as an aid in the diagnosis and should not be interpreted as diagnostics in themselves.

PERFORMANCE CHARACTERISTICS

■ SENSITIVITY

A₄₅₀ of ECP Standard 1 (0 ng/mL) shall be \leq 0.2

A₄₅₀ of ECP Standard 6 (40 ng/mL) shall be \geq 1.2

■ SPECIFICITY

When measuring control samples of known concentrations, each value shall be within $\pm 20\%$ of known concentration.

■ REPRODUCIBILITY

When measuring 3 control samples of known concentrations eightfold simultaneously, CV% of each sample shall be within 15%.

■ INTERFERING SUBSTANCES

Bilirubin F (up to 18.3 mg/dL), Bilirubin C (up to 19.0 mg/dL), Chyle (up to 1,390 unit as Formazine) and /or Rheumatoid factor (up to 500 IU/mL) are not affective on the assay result, but avoid using highly lipemic samples.

Hemoglobin does affect the assay results. Do not use hemolysed samples.

No effect was found to assay values in adding up to 500 ng/mL of EDN.

■ ASSAY RANGE

The assay range of this kit is from 0.125 ng/mL to 40 ng/mL.

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AB Yetkili Temsilcisi: