AllergenA Total IgE

96

RA1002

Reagenti / Reagents / Reagenzien	Quantita/ Quantity/ Menge	
MP	1 x 96	Pronta per l'uso / Ready for use / Gebrauchsfertig
0-5 CAL	8 x 0.5 mL	Pronti per l'uso/ Ready for use / Gebrauchsfertig
CONJ	1 x 18 mL	Pronto per l'uso/ Ready for use / Gebrauchsfertig
CONTROL	2 x 0.5 mL	Pronti per l'uso Ready for use / Gebrauchsfertig
DIL	1 x 3 mL	Pronto per l'uso/ Ready for use / Gebrauchsfertig
SUBS	1 x 14 mL	Pronto per l'uso/ Ready for use / Gebrauchsfertig
WASH 20X	1 x 50 mL	Concentrata 20X / 20X Concentrated / 20x konzentriert
STOP	1 x 14 mL	Pronto per l'uso/ Ready for use / Gebrauchsfertig

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14. LEGENDA SIMBOLI: vedere pagina 39

ENZYME IMMUNOASSAY FOR QUANTITATIVE DETERMINATION OF TOTAL IMMUNOGLOBULIN E (IgE) IN HUMAN SERUM

FOR IN VITRO DIAGNOSTIC USE ONLY FOR PROFESSIONAL USE ONLY

1. CLINICAL APPLICATIONS

Normally IgE concentration in serum is very low. It increases gradually from birth to teen-age. In adults normal concentration of IgE may reach 100 IU/mL. In elderly people IgE level sometimes decreases.

IgE production is essential in anti-helminthic immunity. 15-20-fold increase in IgE concentration is observed in the case of ascariasis. But in industrialized countries detection of high IgE concentrations is mainly connected with allergic diseases. Quantitative determination of total IgE has a great prognostic value. In 75% of children born from parents with allergic diseases, serum IgE concentration is >95% from superior limit of normal range for corresponding age group. Detection of high IgE concentrations in serum by enzyme immunoassay is an important tool for differentiation between allergic diseases and other pathologies with similar clinical manifestations (such as asthma, frequent respiratory diseases, chronic rhinites and dermatites). Increased concentration of total IgE in serum was also reported in patients with lymphosarcoma and Hyper-IgE syndrome.

2. PRINCIPLE OF TEST

AllergenA Total IgE kit is a "sandwich" type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of IgE molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. IgE molecules from the serum sample bind to both immobilized antibody and anti-IgE-peroxidase conjugate (Fig. 1). Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to the IgE level in sample. During incubation with substrate the colour is developing. The intensity of the colour is directly proportional to the concentration of IgE in specimens or calibrators. The IgE concentration in the patient sample is read from a standard curve that is processed in each assay.


Fig. 1. Reaction scheme

3. CONTENT OF THE KIT

MP	Microplate: 1 microplate with 12 breakable 8-well strips (total 96 wells) coated with anti-IgE monoclonal antibodies.
CONJ	Conjugate: 1 vial with 18 mL of solution containing anti-IgE monoclonal antibodies conjugated with HRP. Ready to use
0-5 CAL	IgE Calibrators: 8 vials with 0.5 mL each of protein-based solution containing known IgE concentrations: 0 - 10 - 50 - 100 - 250 - 500 IU/mL . For exact IgE concentrations see vial labels. The set consists of 2 vials of CAL 1 , 2 of CAL 4 and 1 vial respectively: CAL 0 , CAL 2 , CAL 3 and CAL 5 . Ready to use.
CONTROL	IgE Control: 2 vials with 0.5 mL of protein-based solution containing known IgE concentration. For exact range of IgE concentration, see vial labels. Ready to use.
DIL	Sample Diluent: 1 vial with 3 mL of protein-based solution. Ready to use.
WASH 20X	Washing Solution: 1 vial with 50 mL of surfactant in buffered saline, sufficient for preparation of 1000 mL of solution. 20X concentrated
STOP	Blocking Reagent: 1 vial with 14 mL of 1 N HCl solution. Ready to use.
SUB	TMB Solution: 1 vial with 14 mL of 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide. Ready to use.

Note: Extra vials of **CAL 1**, **4** and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for automatic analyser.

3.1 Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Keep **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL **CONTROL** **Calibrators and Control**

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure than no dry matter is left on the caps and walls of the vials.

WASH **20X** Prepare required volume of **Washing Solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH** **20X** + 95 mL of water.

Mix thoroughly, avoid foaming.

SUBS Protect **TMB Solution** from direct light.

3.2 Sample Preparation

Allow samples to reach room temperature (+18...+25 °C). Stir samples gently in order to ensure homogeneity.

If expected IgE concentration in the samples is higher than in calibrator 5, the samples should be diluted **20-fold** with **DIL** **Sample Diluent** in concordance with instructions for use. The example of manual sample dilution as follows:

190 µL of **DIL** **Sample Diluent** + **10 µL** of serum sample,
vortex or mix thoroughly.

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; expiration date for each component is printed on the respective label.

AllergenA Total IgE should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 12 months.

After initial opening the kit is stable for 1 month if stored at +2...+8 °C. If used in several separate experiments, kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with Calibrators, Control, Conjugate, concentrated Washing Solution and Blocking Reagent: at +2...+8 °C until the expiration date;
- vial with TMB Solution: at +2...+8 °C for until the expiration date, protected from light;
- Washing Solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator/shaker (+37 °C, shaking speed 400-800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash)..

6. WARNINGS AND PRECAUTIONS

In order to obtain reproducible results, the following rules must be observed:

- This kit is for *in vitro* diagnostic use only. Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only.

- Do not mix or use together reagents from different lots of the kit except TMB Solution, Blocking Reagent and Washing Solution.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Do not use TMB Solution, Blocking Reagent and Washing Solution supplied by other vendors.
- Use thoroughly clean glassware, free from contamination of metal ions or oxidating substances.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples and solutions; for this purpose, disposable tips should be used for each sample and reagent.
- Do not touch the bottom of the wells.
- **Calibrators** should be measured in each separate assay. It is also recommended to measure each time analyte concentration in the control.
- Reconstitute lyophilized reagents, if present, as described on the relative labels and Instructions for use. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
- If kit is used in several separate runs, it is necessary to take reagents from analyzer immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator.
- Recalibration using calibration curve, obtained with kit of any other lot, is not allowed.
- TMB Solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
- Utilize a suitable method for the correct identification of patient samples. Incorrect identification may cause a specificity losses of the system and wrong clinical results.

In order to avoid personal and environmental contamination, the following precautions must be observed:

- Use protective individual devices (ex.: disposable gloves, lab coats, etc.) while handling potentially infectious material as well as during the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- TMB Solution and Blocking Reagent (1 N HCL) should be handled with care, as well as the reagents containing Kathon (Conjugate, Calibrators, Control and Sample Diluent). Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.

- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious. Therefore, the assay waste must be disposed of in accordance with established safety procedures.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.

As the kit contains irritant (**CONJ E** **DIL** **CAL** **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

- According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code 180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

7. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipemic (milky) serum samples as well as samples containing sodium azide as preservative. Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze samples for longer storage (-20 °C and lower). Avoid repeated freezing.

8. TEST PROCEDURE

The manual procedure must be performed as described in paragraph 8.2. For internal quality control, it is advisable to use Control Serum at known concentration that allow assay performances to be monitored.

8.1 Assay Procedure for Automatic Test

We guarantee its applications on RADIM and/or SEAC, Next Level automatic instruments.

While using another automatic instrument for microplate, it is under end user's responsibility to make sure that it was appropriately tested for ELISA kits.

While using for the procedure automatic analyser, refer to its relative manual.

Note: *If the kit is used on instrument in several separate experiments, it is necessary to take reagents out of analyzer immediately after pipetting them into the wells of all plates, because liquid evaporates from vials. Put the reagents into refrigerator.*

8.2 Assay Procedure for Manual Test (see assay scheme, page 40)

A. Pipette: **150 µL** of **Conjugate** **CONJ** into each well **except wells A1-A2 (blank)**.

B. Pipette: **20 µL** of **Calibrators** **CAL**, **Control** **CONTROL** and patient's samples in duplicates into the respective wells.

Leave wells A1-A2 empty for blank.

C Incubate strips for **90 minutes at +37 °C while shaking (500-800 rpm)**.

D. Wash **5 times**, as described in section 8.3.

E. Pipette **100 µL** of **TMB Solution** **SUBS** into each well (including blank); incubate strips **at room temperature (+18...+25 °C) in the dark for 15-30 minutes**, depending on the colour intensity, or **10 minutes while shaking (500-800 rpm) at +37 °C**.

F. Pipette **100 µL** of **Blocking Reagent** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing **TMB Solution**. Shake for **1-2 min at room temperature**.

G. Read OD at **450 nm within 20 min**.

8.3 Wash Procedure

It is advisable to use an automatic microplate washer set at **5 wash cycles** and a volume of **300 µL** of **wash solution** per well per cycle. If an automatic washer is not available, the wash procedure can be carried out manually:

- remove the contents of the wells into a container with disinfectant;
- dispense **300 uL** of wash solution (prepared according to section 3.1) into each well, shake the plate carefully for **5-10 sec** and remove the contents of the wells; repeat **5 times**;
- strike the wells sharply on absorbent material to remove any liquid residue.

9. DATA PROCESSING

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = $2.28 - 0.06 = 2.22$

9.1 Data Reliability (for OD Measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) <0.100;
- average OD of Cal 5 >1.0 (after blank subtraction).
- control's concentration must fall within the acceptability range that is shown on the vial label.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

9.2 Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective IgE concentrations using 4PL or 5PL fit (see typical standard curve, fig. 2). Calculate concentration of IgE in samples using standard curve.

Any extrapolation of the standard curve to IgE concentration above the nominal value of the calibrator 5 (approximately 500 IU/mL) is forbidden. In this case the sample should be diluted 20-fold with **Sample Diluent** and re-tested. Multiply the measured concentration of pre-diluted samples by dilution factor (20-fold).

If concentration of the sample diluted 20-fold remains above concentration of the **CAL** 1-5 to assign concentration "above 10000 IU/mL" to a sample.


Fig. 2. Example of typical standard curve.

Do not use for evaluation of real assay data.

10. EXPECTED VALUES

Serum samples collected between 9 and 11 a.m. from 79 apparently healthy people (both males and females) between the ages of 21 -45 and 226 patients with allergic diseases (of the same age), were assayed with **AllergenA Total IgE** kit. Among healthy individuals in 48% cases IgE concentration was below 25 IU/mL; in 34% cases between 25-100 IU/mL, and only in 18% cases above 100 IU/mL. Among allergic patients IgE concentration below 25 IU/mL was detected in 8% cases; between 25-100 IU/mL in 22% cases and above 100 IU/mL in 70% cases. These limits should be considered as guidelines only.

It is highly recommended for each laboratory to determine its own reference range of IgE concentrations.

11. QUALITY CONTROL

The use of control samples is advised to assure the day to day validity of results.

12. PERFORMANCE CHARACTERISTICS OF THE ASSAY

12.1 Calibration-Traceability

AllergenA Total IgE kit was calibrated against the WHO 2nd International Reference Preparation 75/502.

12.2 Specificity

According to Product specification, provided by the supplier, no cross-reactivity was detected between anti-IgE monoclonal antibodies used in the assay and IgG, IgM and IgA.

12.3 Analytical Sensitivity

Analytical sensitivity of **AllergenA Total IgE** kit, i.e. concentration that can be distinguished from zero calibrator, is 2.3 IU/mL. It is defined as mean OD of 10 replicates of Calibrator 0 plus 2 SD.

12.4 Measurement Range

AllergenA Total IgE kit was validated for measurement of IgE concentration within the concentration range (without dilution) of 2.3 - 500 IU/mL.

12.5 Measurement Units

In **AllergenA Total IgE** kit the concentrations of calibrators are specified in IU/mL. To convert into ng/mL, multiply the concentration in IU/mL by 2.4.

12.6 Hook Effect

For **AllergenA Total IgE** kit **high dose hook effect** was not detected for concentrations up to 10 000 IU/mL. **High dose hook effect** was determined by spiking Calibrator 0 with antigen.

12.7 Intra- and Inter-Assay Variation (Precision)

For **intra-assay CV** determination, 8 serum samples were run, each in 9 replicates. The results are showed below.

Sample	Mean IgE concentration, IU/mL	Intra-assay CV	
		SD	CV, %
HS 1	14.1	1.01	7.2
HS 2	38.5	1.64	4.3
HS 3	81.4	4.35	5.4
HS 4	88.0	5.54	6.3
HS 5	103	4.7	4.5
HS 6	195	5.8	3.0
HS 7	262	18.1	6.9
HS 8	434	22.9	5.3

For **inter-assay CV** determination 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are showed below.

Sample	Mean IgE concentration, IU/mL			Inter-assay CV	
	Assay 1	Assay 2	Assay 3	SD	CV, %
HS 1	14.1	15.0	14.6	0.45	3.1
HS 2	38.5	37.5	41.2	1.91	4.9
HS 3	81.4	82.1	82.7	0.65	0.8
HS 4	88.0	84.6	89.4	2.47	2.8
HS 5	103	99	106	3.8	3.7
HS 6	195	185	208	11.4	5.8
HS 7	262	239	272	17.0	6.6
HS 8	434	398	423	18.6	4.4

13. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings. Due to the heterogeneous nature of human antibodies there might be samples that do not maintain dilution parallelism.

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